



## Cleanse (no rinse)

For frequent cleansing and moisturising when soap and water may otherwise dry the skin

# TENA Wet Wipes

These 3-in-1 soft and thick wipes are pre-moistened to cleanse, restore and protect in one simple step.

- Primarily used for frequent perineal skin care of incontinent residents when changing absorbent products, but can also be used for full body cleansing
- One-handed dispensing and 'clean touch' system avoids the risk of touching the next wipe in the pack to minimise the risk of cross contamination
- A single wipe can cleanse a large area of skin before disposal
- Helps to maintain the skin's natural pH, dermatologically tested, and mildly perfumed



Product Name	Order #	Description	Pcs/Case	Format
TENA Wet Wipe	6479	-	576	12/48
TENA Wet Wipe	6511	-	600	60/10

Ingredient (INCI name)	Role
Aqua	Solvent/carrier
PEG-7 Cocoglycerides	Emollient/ Skin conditioning/ cleansing agent
PEG-40 Hydrogenated Castor Oil	Emollient/ Skin conditioning/ cleansing agent
Phenoxyethanol	Preservative
Tetrasodium EDTA	Chelating agent
Benzoic Acid	Preservative
Glycerin	Emollient/ Skin conditioning/ cleansing agent
Glyceryl Oleate	Emollient/ Skin conditioning/ cleansing agent
Dehydroacetic Acid	Preservative
Polyaminopropyl Biguanide	Preservative
Ethylexylglycerin	Preservative
Parfum	Perfuming

## Safety

The product fulfils the essential safety requirements of Regulation (EC) No. 1223/2009 of the European Parliament and the Council of 30 November 2009 on cosmetic products. As stipulated in Article 18 the product has not been tested on animals. The production unit is following the international standard ISO 22716:2007 as regards product documentation and Good Manufacturing Practices, GMP. The product has been evaluated with epicutaneous testing under dermatological control.

## First aid

Unless an excessive amount is swallowed, the impregnating fluid is unlikely to be harmful. If ingestion does occur, remove material from the mouth and drink water. If discomfort is experienced after eye contact, rinse with water. Consult a physician in case of problems due to prolonged and heavy exposure.

## Physical properties

The wipe material is an embossed spunlaced material of nonwoven of polyester and viscose fibres impregnated with a lightly scented emulsion

**Size:** 29 cm x 22 cm

**pH:** 5.5

## Storage and period of use

### Product durability after production

**(expiry date):** 30 months

### Recommended storage temperature:

Ambient temperature between 5°C and 25°C

## For further information, please contact:

hcmarketing@sca.com

SCA Hygiene Products AB  
[www.tena.com](http://www.tena.com)



**TENA Skincare**

*Together we make a difference*